

## **PerCP Conjugation Kit**

Cat No: BIOK1968

Description

PerCP (Pyrophyllin-Chlorophyll-Protein Complex) is a new type of fluorescent labeling reagent commonly used at present, with a maximum excitation wavelength of 472nm and a maximum emission wavelength of 677nm. It is used for multi-label flow technology and biochip detection. It has stable properties: Fluorescence quantum yield is high, background interference is small, and it is easy to cross-link with biotin antibodies and glycoproteins. As a new generation of fluorescent probes, it can replace isotopes and enzymes in clinical diagnosis, immunology, molecular biology as a marker.

Biotyscience's PerCP antibody labeling kit is used for rapid labeling of PerCP and antibodies. The labeling can be completed within 4 hours. It is easy to operate, at least 100ug of antibody can be labeled, 100% of the sample is recovered, and the labeled product is stable between batches and has good repeatability.

## **PRODUCT INFORMATION**

Components	100ug	500ug	1mg	Preservation conditions
Lightning-Link label (Lyophilized)	1 vial	1 vial	1 vial	-20 $^{\circ}$ , Protect from light
LL-Modifier reagent	100ul	100ul	150ul	-20 $^{\circ}$ C, Protect from light
LL-Quencher reagent	100ul	100ul	150ul	-20 $^{\circ}$ C, Protect from light

Antibody is in liquid state, dissolved in 1×PBS, concentration 1mg/ml.

北京百欧泰生物科技有限公司 Tel: +86 136 8125 6816 Email: biotyscience@163.com Address: 北京市房山区良乡凯旋大街建设路 18 号-D3032



Note: If necessary, the antibody needs to be dialyzed again or passed through a desalting column for accurate protein quantification.

## Steps

1. According to the ratio: add 100ul antibody (100ug) to 10ul LL-Modifier and mix well.

2. Add the mixed antibody to the corresponding labeled amount of Lightning-Link label dry powder, and react at room temperature for 3 hours or overnight in the dark.

3. According to the ratio: add 100ul antibody (100ug) to 10ulLL-Quencher, mix well, and complete the labeling after reaction for 30min.

4. Store at 2-8°C in the dark, do not freeze. The antibody can be used directly without adding other ingredients.

## **Contact Us**

Beijing Biotyscience Co. Ltd.

**QQ**: 499854788

3494243873

1619694293

We chat: 13681256816

15511114213

Email: biotyscience@163.com

Tel:+86 136 8125 6816