# **R-Phycoerythrin Antibody Conjugation Kit**

#### Cat No: BIOK1969

Description

Phycoerythrin (R-Phycoerythrin, R-PE) is separated and purified from red algae. It is a new type of fluorescent labeling reagent commonly used at present. The excitation wavelength is 564nm and the emission wavelength is 578nm. It is used for multi-label flow technology and biochip detection. It has the characteristics of stable properties, high fluorescence quantum yield, low background interference, and easy cross-linking with macromolecules such as biotin antibodies and glycoproteins. The new generation of fluorescent probes can replace isotopes and enzymes as markers in clinical diagnosis, immunology, molecular biology, etc.

Biotyscience' R-PE antibody labeling kit is used for rapid labeling of R-PE and antibodies. The labeling can be completed within 4 hours. The operation is simple. At least 100ug of antibody can be labeled. 100% of the sample is recovered. The labeled product is stable between batches. Good repeatability.

### **PRODUCT INFORMATION**

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Components	500ug	1mg	Preservation conditions
Lightning-Link label (Lyophilized)	1 vial	1 vial	-20℃, Protect from light
LL-Modifier reagent	100ul	150ul	-20°C, Protect from light
LL-Quencher reagent	100ul	150ul	-20 <sup>°</sup> C, Protect from light

### Steps

#### 1. Antibody preparation

1.1.Antibody is in liquid state, dissolved in  $1 \times PBS$ , concentration 1mg/ml Note: If necessary, the antibody needs to be dialyzed again or passed through a desalting column for accurate protein quantification

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**Beijing Biotyscience Co. Ltd.**