
FITC-Antibody Conjugation Kit

Cat No: B1OK1967

Description

Fluorescein-labeled antibody is abbreviated as fluorescent antibody (FA), which is currently widely used in immunopathology, cytochemistry, flow cytometry, virology and autoantibody clinical immunodiagnosis of specific, sensitive, qualitative and localized immune Chemical reagents. Antibodies used for labeling require high specificity and high affinity. The antiserum used should not contain antibodies against normal tissues in the specimen. Generally, IgG and IgM need to be purified and extracted before labeling.


The fluorescein in this kit uses high-quality imported fluorescein isothiocyanate (FITC). The principle of FITC-labeled antibody is to use the $\text{N}=\text{C}=\text{S}$ group on FITC and the free $-\text{NH}_2$ on the antibody to chemically form FITC-Antibody conjugates. There are 86 lysine residues in an IgG molecule, and generally can bind 15-20 at most, and an IgG molecule can bind 2-8 molecules of FITC.

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NOTES

1. Store the kit at 2-8℃, do not freeze.
2. The components of the kit may be upside down during transportation, causing liquid or dry powder reagents to stick to the tube wall or bottle cap. Please centrifuge before use to make the liquid or dry powder reagent attached to the tube wall or bottle cap settle to the bottom of the tube.
3. FITC needs to be prepared for immediate use, and the dissolved FITC cannot be stored for a long time.
4. The buffer in the purification column contains the toxic component sodium azide

(NaN₃). Avoid contact with skin, eyes and mucous membranes when using it.

5. DMSO is slightly toxic, permeable to human skin and irritating to the eyes. Avoid contact with skin, eyes and mucous membranes when using it.

6. The dialysis, concentration and concentration determination of the antibody before labeling will all cause the loss of the amount of antibody. Therefore, when preparing the antibody before labeling, consider the most appropriate amount of antibody according to the specific situation.

Contact Us

Beijing Biotyscience Co. Ltd.

QQ: 499854788

3494243873

We chat: 13681256816

15511114213

Email: info@biotyscience.com

Tel: 400-669-8850