Cy3-Antibody Conjugation Kit

Cat No: BIOK1963

Description

Immunofluorescence technique is to label a known antibody or antigen molecule with fluorescein. When it reacts with the corresponding antigen or antibody, there is a certain amount of fluorescein on the formed complex. Below you can see the fluorescent antigen-antibody binding site and detect the antigen or antibody. The main features of this technology are: strong specificity, high sensitivity and fast speed.

Fluorescein-labeled antibody is abbreviated as fluorescent antibody (FA), which is currently widely used in immunopathology, cytochemistry, flow cytometry, virology and autoantibody clinical immunodiagnosis of specific, sensitive, qualitative and localized immune Chemical reagents. Antibodies used for labeling require high specificity and high affinity. The antiserum used should not contain antibodies against normal tissues in the specimen. Generally, IgG and IgM need to be purified and extracted before labeling.

The fluorescein in this kit uses high-quality imported Cy3. The cyanine dye Cy3, SE is one of the few near-infrared fluorescent labeling dyes. It produces a low red fluorescent background, high fluorescence intensity, and good fluorescence stability. Cy3, SE are often used in biomolecule labeling, fluorescence imaging and other fluorescent biological analysis, such as labeling peptides, proteins and nucleic acids. Cy3, SE fluorescent labeling dye can significantly enhance the fluorescence intensity after binding to protein.

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
NOTES

1. Store the kit at 2-8℃, do not freeze.
2. The components of the kit may be upside down during transportation, causing liquid or dry powder reagents to stick to the tube wall or bottle cap. Please centrifuge before use to make the liquid or dry powder reagent attached to the tube wall or bottle cap settle to the bottom of the tube.
3. The buffer in the purification column contains the toxic component sodium azide (NaN₃). Avoid contact with skin, eyes and mucous membranes when using it.
4. DMSO is slightly toxic, permeable to human skin and irritating to the eyes. Avoid contact with skin, eyes and mucous membranes when using it.
5. The dialysis, concentration and concentration determination of the antibody before labeling will all cause the loss of the amount of antibody. Therefore, when preparing the antibody before labeling, consider the most appropriate amount of antibody according to the specific situation.

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Beijing Biotyscience Co. Ltd.